
FINGER INJURY SYSTEM

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ACU-LIFE®

Every injury is different and requires specific treatment. Please consult your physician prior to starting treatment.

Suggestions for treating a jammed finger:

1. Use the gel pack for cold therapy 3-4 times a day for 10 minutes until the swelling and pain subsides. Slide gauze tube over finger before applying cold pack. Place cold pack around finger and use the strap to fasten the cold pack to finger (Note Photo A). Strap should not be positioned where it will put direct pressure on the knuckle or any wound.

2. Wear the comfort splint to help protect the finger while it heals. Use the strap to fasten the splint to the finger (Note Photo B). Splint should feel snug but not too tight. Strap should not be positioned where it will put pressure directly on the knuckle or any wound. Strap can be removed and repositioned if necessary. If redness, pain or swelling persists, consult a physician at once.

Finger Cold Pack Use:

Place gel pack in freezer for at least 2 hours prior to using. Always return to freezer promptly for fastest reuse.

Finger Warm Pack Use:

Place gel pack in warm water for 10 minutes. Do not use the microwave or boiling water as this will over-heat the gel pack and cause serious injury.

Cleaning:

Wash gel pack in lukewarm water and lay flat to dry.

Warning:

Cold therapy should not be used by individuals with known or suspected circulation problems. Do not apply heat to broken or sensitive skin. Heat should not be used with external analgesic products. Do not heat gel pack in a microwave. Do not sit or lean on gel pack, excessive weight may cause the pack to rupture and leak.

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Health Enterprises, Inc.