

ACU-LIFE®




THERAPEUTIC FOOT MASSAGER

Item# 400708

Helps Relieve Pain & Reduce Stress

- A** *Therapeutic Massage Knobs* allow user to control amount of pressure for a gentle or deep tissue massage.
- B** *No-Slip Rings* provide traction for maximum control on hard surfaces (hardwood, tile, marble, etc.).
- C** *Raised Edge* design secures foot to help prevent slippage during use.
- D** *Removable Top* makes it easy to add ice and cold or warm water for a targeted hot or cold massage.



Instructions:

- Sit down and place Massager in the center of the arch.
- While sitting with your knee bent to 90°, gently roll your foot over the massager.
- Apply pressure onto massager and roll back and forth.

Cold Therapy:

- Unscrew cap from one end of massager.
- Add ice cubes and cold water.
- Replace cap and begin massage.

Caution: Do not place massager in freezer. Make certain that cap is tightly secured to prevent leakage.

Warm Therapy:

- Unscrew cap from one end of massager.
- Add warm water.
- Replace cap and begin massage.

Caution: Do not place massager in microwave or use boiling water. Be careful when filling massager to prevent skin contact with warm water. Make certain that cap is tightly secured to prevent leakage.

Caution:

If you have any serious medical conditions, consult your physician prior to use. If you experience pain during use, decrease the amount of pressure or treatment time. If symptoms continue, discontinue use and consult your physician.

Cleaning Instructions:

Clean the Hot & Cold Foot Massager with warm soapy water and wipe with a clean dry cloth or air dry. After each use remove top (D), carefully empty all water and dry inside of Massager.

HE® Health Enterprises, Inc.
90 George Leven Drive
North Attleboro, MA 02760
www.healthenterprises.com
info@healthenterprises.com

For more unique healthcare products, visit us at www.healthenterprises.com

Also Try...

ACU-LIFE®

24/7 PLANTAR FASCIITIS

Rev. 072208a

Item# 400706



Night Splint

Stretches Plantar Fascia & Achilles Tendon



Day Arch Brace
Increases Comfort While Standing & Walking

TRU-ICE™

Full Benefit of an Ice Bag or Cold Compress in 1/3 the Time



Plastic Liner Prevents Mess

Ice Bag or Cold Compress
15-20 minutes

Tru-Ice™
5-7 minutes

For Carpal Tunnel, Wrist Sprains, Strains, Arthritis & Tendonitis:



For Ankles, Elbows and Knees:



For Stiff, Strained or Aching Backs:



HE® Health Enterprises, Inc.
90 George Leven Drive
North Attleboro, MA 02760
www.healthenterprises.com
info@healthenterprises.com

For more unique healthcare products, visit us at www.healthenterprises.com