

75% of all hearing aid repairs are due to moisture and ear wax.

Almost 100% of these repairs are preventable with proper daily maintenance.

(Source: healthyhearing.com)

Earwax Problem

- Hearing Aids can block earwax from exiting the ear canal and lead to earwax build-up inside the hearing aid.
- Earwax is a very dense material that can impact sound quality and lead to costly repairs.

Moisture Problem

- The typical ear canal is a warm, humid environment, which creates moisture that can damage the sensitive electronics of a hearing aid.
- Unlike earwax, moisture in a hearing aid is difficult to see because it is typically caused by moisture molecules, not visible droplets.

Earwax & Moisture Can Cause:

- Whistling or feedback.
- Intermittent (on/off) or distorted sound.
- Costly and inconvenient repairs.
- Reduced hearing aid life.

Hearing Aid Maintenance Tips

- Always make sure that your hands are dry and clean before handling your hearing aid.
- Always clean your hearing aid from below to allow wax or dirt to fall out instead of going back inside the hearing aid.
- Clean your hearing aid over a soft cloth on top of a table or desk to reduce the risk of damage if the hearing aid is accidentally dropped.
- Be careful to follow the cleaning instructions provided by your hearing healthcare professional.
- When you are not using your hearing aid, turn it off, remove the battery and place the hearing aid into the Dri-Eze™ Dehumidifier.
- Never insert any sharp tools into the microphone opening as damage can occur.



HEARING AID CLEANING WIPES

Item# 400675

- Removes wax and debris from hearing aids and earmolds.
- 30 individually wrapped pre-moistened wipes.



HEARING AID BATTERY TESTER

Item# 400692

- Quickly and easily determine if your hearing aid batteries are still good.
- Works on all hearing aid batteries.



HEARING AID AIR BLOWER

Item# 400674

- Helps keep earmolds and tubing clear from moisture, dirt and debris.
- Improves sound quality.



HEARING AID PHONE EARPADS

Item# 400676

- Reduces feedback and whistle noise.
- Improves listening comfort.
- Contains two universal phone pads and one BONUS cell phone pad.

ACU-LIFE®

DAILY HEARING AID CARE & MAINTENANCE:



- Removes Damaging Earwax & Moisture
- Improves Sound Quality
- Prevents Costly Repairs
- Prolongs Hearing Aid Life
- Works With All Hearing Aids



Health Enterprises, Inc.
90 George Leven Drive
North Attleboro, MA 02760
info@healthenterprises.com
www.healthenterprises.com

Call with questions and comments:
1.800.633.4243